



Laboratory rats, mice and rodent capture for biomedical research

- Multiple camera setup.
- Compatible with GigE Vision, USB3, 10 and 25 GigE, Camera Link and CoaXPress cameras.
- Supports cameras and grabbers from over 200 different OEM's using native API's.
- Synchronize camera exposure over multiple cameras and over multiple computers.

Capture and synchronize video with DAQ sensor data such as :

- ECG, EEG, Temperature, etc...

Supports cameras with IR and NearIR Light sensors :

- Capture in Pre-Post using triggering to save disk space.
- Monitor over URL using Remote.



NORPIX

DIGITAL VIDEO RECORDING SOFTWARE

StreamPix - TroublePix - StreamTouch

► www.norpix.com



Laboratory rats, mice and rodent capture for biomedical research

- Multiple camera setup.
- Compatible with GigE Vision, USB3, 10 and 25 GigE, Camera Link and CoaXPress cameras.
- Supports cameras and grabbers from over 200 different OEM's using native API's.
- Synchronize camera exposure over multiple cameras and over multiple computers.

Capture and synchronize video with DAQ sensor data such as :

- ECG, EEG, Temperature, etc...

Supports cameras with IR and NearIR Light sensors :

- Capture in Pre-Post using triggering to save disk space.
- Monitor over URL using Remote.



NORPIX

DIGITAL VIDEO RECORDING SOFTWARE

StreamPix - TroublePix - StreamTouch

► www.norpix.com